



Abraxis LLC

Malachite Green Plate Kit

PN 52251B

Instructional Booklet

Read Completely Before Use.

INTENDED USE

The Abraxis Malachite Green plate kit is a competitive ELISA for the quantitative analysis of Malachite green in aquatic products.

ASSAY PRINCIPLES

The Abraxis Malachite Green ELISA plate kit is a competitive enzyme-labeled immunoassay. Malachite Green is extracted from a sample by blending or shaking with extraction solvent. The diluted Malachite Green sample extract and calibrators are pipetted into test wells followed by malachite green antibody to initiate the reaction. After a 30 minute incubation, the wells are washed with laboratory grade water. Malachite Green HRP conjugate is then added and incubated for 30 minutes. Following a second wash after the incubation, a clear substrate is added to the wells and any bound enzyme conjugate causes the conversion of the colorless substrate/chromogen to a blue color. After a 20-30 minute incubation, the reaction is stopped and the amount of color in each well is read. The color of the unknown samples is compared to the color of the calibrators and the Malachite green concentration of the samples is determined by interpolation.

SPECIFICITY

The Abraxis Malachite green Plate Kit can not differentiate between the various Malachite Greens metabolites, but detects their presence to differing degrees. The following table shows the % cross reactivity of Leucomalachite versus Malachite green.

Compound	%CR
Malachite green	100%
Leucomalachite green	<1%
Crystal Violet	120%
LeucoCrystal Violet	<1%

DETECTION LIMIT

Aquatic products: 0.25ppb



Abraxis LLC

REAGENTS AND MATERIALS PROVIDED

The kit in its original packaging can be used until the end of the month indicated on the box label when stored at 2-8°C.

- Plate containing 12 test strips of 8 wells each vacuum-packed in aluminized pouch with indicating dessicant.
- 1 vials containing 0.4 ml of Malachite Green calibrator stock solution 100µg/L(ppb) of Malachite Green.
- 1 vial containing 0.8 ml Malachite Green HRP Enzyme Conjugate concentrate X 25.
- 1 vial containing 9 ml of anti-Malachite Green antibody.
- 1 vial containing 16 ml of Substrate.
- 1 vial containing 12 ml of Stop Solution (Caution 1N HCl . Handle with care.)
- 1 vial containing 15ml HRP conjugate diluent.
- 2 X 25 mL bottles of Sample Dilution Buffer
- 1 bottle containing 100 mL of Wash Solution concentrate 5X. Needs dilution prior to use.
- Reagent A (4X) (for Sample preparation)15 ml. **Remove from kit, aliquot and freeze upon receipt of kit.**
- Reagent B (for Sample preparation)12 g.
- Alumina , Neutral 100-200 mesh, 130 g.
- Instructions.

PRECAUTIONS

1. Each reagent is optimized for use in the Abraxis Malachite Green Plate Kit .Do not substitute reagents from any other manufacturer into the test kit. Do not combine reagents from other Abraxis Malachite Green Plate Kits with different lot numbers.
2. Dilution or adulteration of reagents or samples not called for in the procedure may result in inaccurate results.
3. Do not use reagents after expiration date.
4. Reagents should be brought to room temperature, 20-28°C (62-82°F) prior to use. Avoid prolonged (>24hours) storage at room temperature.
5. Malachite green is an antibiotic and should be treated with care.
6. For longer time storage, Reagent A should be removed from kit upon arrival, aliquoted (2 mL/vial) and frozen at – 20 °C. One of the frozen vials can be used for dilution (**Preparation of Reagents section**).
7. The Stop Solution is 1N hydrochloric acid, avoid contact with skin and mucous membranes. Immediately clean up any spills and wash area with copious amounts of water. If contact should occur, immediately flush with copious amounts of water.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Laboratory quality distilled or deionized water.
2. Graduated cylinder,100 ml or larger .
3. Glassware for sample extraction and extract collection.
4. Acetonitrile.



Abraxis LLC

5. Pipet with disposable tips capable of dispensing 50 μ L
6. Multi-channel pipet; 8 channel capable of dispensing 50 and 100 μ L
7. paper towels or equivalent absorbent material
8. Microwell plate or strip reader with 450nm filter
9. Timer.
10. Vortex mixer.

Preparation Reagents

Reagent A: Dilute Reagent a 1:4 with Acetonitrile (100%). Prepare fresh before each assay.

Standard dilution buffer: Mix 1 mL of acetonitrile with 9 mL of sample diluent provided. Cover and swirl to mix completely

HRP-conjugate solution: The conjugate provided in the kit is 25X concentrate. Before each assay, the conjugate should be diluted into conjugate diluent. For example mix 100 μ L of 25X concentrate conjugate with 2.4 mL of conjugate diluent to obtain the working conjugate.

CALIBRATOR PRAPARATION

The Malachite green calibrator stock solution is diluted in **Standard dilution buffer** prior to use in the assay. Bring the stock Malachite Green calibrator stock solution to room temperature prior to diluting.

The stock standard will need to be diluted daily in 10% acetonitrile/sample dilution buffer (v/v) to prepare calibrators corresponding to 0, 0.005, 0.01, 0.02, 0.05, 0.2, and 0.5 μ g/L (ppb) of Malachite Green.

Label 6 clear 4 mL glass tubes or vials: “0.5”; “0.2”; “0.05”; “0.02”; “0.01”; “0.005”

- Standard at 0.5 ppb: 20 μ L of 100 ppb stock into 3.98 mL of standard diluent
- Standard at 0.2 ppb: 400 μ L of the 0.5 ppb standard into 600 μ L of standard diluent
- Standard at 0.05 ppb: 200 μ L of the 0.5 ppb standard into 1.8 mL of standard diluent
- Standard at 0.02 ppb: 400 μ L of the 0.05 ppb standard into 600 μ L of standard diluent
- Standard at 0.01 ppb: 400 μ L of the 0.05 ppb standard into 1.6 mL of standard diluent.
- Standard at 0.005 ppb: 500 μ L of the 0.01 ppb standard into 500 μ L of standard diluent

Be sure to tightly cap the vials containing the stocks to minimize evaporative losses.

SAMPLE PREPARATION (Meat 1:50 Dilution)

1. Weigh 5g of tissue into a 20 or 40 mL glass vial
2. Add 10mL of Acetonitrile
3. Vortex (at max speed) for 1 minute



Abraxis LLC

4. Place vial(s) in shaker/rotator for 10 minutes
5. Centrifuge vial(s) for 10 minutes at 4,000rpm
6. Extract all supernatant into new vial(s)
7. Add 3g of Alumina (neutral) to vial(s)
8. Vortex for 1 minute
9. Centrifuge vial(s) for 10 minutes at 4,000rpm
10. Extract 2mL of supernatant into new vial(s)
11. Evaporate supernatant at 60°C
12. Add 0.5mL of Reagent A (diluted 1:4) to dried vial(s)
13. Vortex for 30 seconds
14. Incubate vial(s) for 30 minutes
15. Add 0.1g of Reagent B to vial(s)
16. Vortex for 10 seconds and let stand for 15 minutes
17. Dilute 100 μ L of supernatant with 1.9 mL of standard diluent (sample diluent buffer with 10% Acetonitrile)
18. Dilute 1:5 in sample diluent
19. Use 75 μ L of sample in assay
20. Multiply the result by 50 to account for dilution.

TEST PROCEDURE (Note: Running calibrators and samples in duplicate will improve assay precision and accuracy.)

1. Allow reagents and sample extracts to reach room temperature prior to running the test. Fill a wash bottle with wash solution.
2. Place the appropriate number of test wells and into a microwell holder. Be sure to re-seal unused wells in the zip-lock bag with dessicant.
3. Using a pipet with disposable tips, add **75 μ L of Calibrators or Sample extract** to the appropriate test wells. Be sure to use a clean pipet tip for each. Add **75 μ L of Antibody solution** to each well.
4. Rotate the plate gently for 30 seconds and incubate the test wells for **30 minutes**.
5. Decant the contents of the wells into an appropriate waste container. Fill the wells to overflowing with wash solution and decant. Repeat 3X for a total of four washes. Tap the plate against absorbent paper.
6. Dispense **150 μ L of Enzyme Conjugate** into each test well.
7. Rotate the plate gently for 30 seconds and incubate the test wells for **30 minutes**.
8. Decant the contents of the wells into an appropriate waste container. Fill the wells to overflowing with wash solution and decant. Repeat 3X for a total of four washes. Tap the plate against absorbent paper.
9. Dispense **150 μ L of Substrate** into each well.
10. Incubate the wells for **20-30 minutes**.
11. Dispense **100 μ L of Stop Solution** into each test well.
12. Read and record the absorbance of the wells at 450nm using a strip or plate reader.



Abraxis LLC

RESULTS INTERPRETATION

1. Semi-quantitative results can be derived by simple comparison of the sample absorbances to the absorbance of the calibrator wells: Samples containing less color than a calibrator well have a concentration of Malachite Green greater than the concentration of the calibrator. Samples containing more color than a calibrator well have a concentration less than the concentration of the calibrator.
2. Quantitative interpretation requires graphing the absorbances of the calibrators (Y axis) versus the log of the calibrator concentration (X axis) on semi-log graph paper. A straight line is drawn through the calibrator points and the sample absorbances are located on the line. The corresponding point on the Y axis is the concentration of the sample. Samples with absorbances greater than the lowest calibrator or less than the highest calibrator must be reported as <0.25 ppb or >25 ppb, respectively.

Alternatively, Abraxis can supply a spreadsheet template which can be used for data reduction. Please contact Abraxis for further details.

GENERAL LIMITED WARRANTY

Abraxis LLC ("Abraxis") warrants the products manufactured by it against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product's printed expiration date. ABRAXIS MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of Abraxis products appearing in published catalogues and product literature may not be altered except by express written agreement signed by an officer of Abraxis. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and, if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Abraxis's sole obligation shall be to repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period; provided the customer notifies Abraxis promptly of any such defect. The exclusive remedy provided herein shall not be deemed to have failed of its essential purpose so long as Abraxis is willing and able to repair any nonconforming Abraxis product or part. Abraxis shall not be liable for consequential, incidental, special or any other indirect damages resulting from economic loss or property damage sustained by a customer from the use of its products. However, in some states the purchaser may have rights under state law in addition to those provided by warranty.

R090911

ABRAXIS LLC
54 Steamwhistle Drive
Warminster, PA 18974

Tel. (215) 357-3911
Fax (215) 357-5232
WEB: www.abraxiskits.com